

AE-TU200 Benchtop Turbidity Meter

Features:

- ◆ 1-5 points calibration using Formazin standards or zero water.
- ◆ Autoranging from 0-200 NTUs.
- ◆ 850nm LED Light Source.
- ◆ Uses Nephelometric principle with 90° detection.
- ◆ Meets ISO 7027 and offers consistent light source for higher accuracy.
- ◆ Selectable turbidity units, including the NTU, FNU and EBC.
- ◆ Settable parameters, including user ID, number of calibration points, data and time, etc.



Measuring slot
Convenient flip design



Meter includes
5pcs turbidity bottles

Technical Data

Model		TU200
Measuring Unit	Range	(0 ~ 20.00) NTU (20.0 ~ 200.0) NTU
	Accuracy	≤ ±8%
	Repeatability	≤1.0%
	Fluctuation	≤ ±1.0% FS/30min
General	Power Supply	AC Adapter, 200-240 VAC input
	Dimension (mm)& Vtetcht (kg)	330×250×130, 2.5