

Disintegration Tester



ADVANTAGE:

- Auto-diagnose , auto-alarm.
- Stainless steel basket and paddle.
- Micro processor and LCD Display.
- Each beaker can work alone.
- Water circulating in the bath to keep the temperature conformity with magnetic pump.

Specification

Model	BJ-I	BJ-II	BJ-III	BJ-IV
Basket	1 pc	2 pcs	3 pcs	4pcs
Speed Range	30~31rpm			
Basket Lifting Height	55±1mm			
Minimum Spacing	25±2mm			
Mesh Size	Standard 2mm (0.425mm, 1mm, 710UM apertures available)	Standard 2mm (0.425mm, 1mm, apertures available)	Standard 2mm (0.425mm, 1mm, 710UM apertures available)	Standard 2mm (0.425mm, 1mm, apertures available)
Temp. Range	RT~50°C			
Temp. Precision	±0.3°C			
Timing Range	1~900min			
Time Precision	±0.5min			
Power Consumption	600W			
Power Supply	AC220V/110V±10%, 50/60HZ			
External Size(W*D*H)	260*300*440 mm	400*320*440 mm	540*320*440 mm	400*320*440 mm