

Microscope XSZ-107BN

Specification

45mm DIN Achromatic Objective:4X,10X,40X(S),100X(S,Oil)

Wide Field Plan-Scope Eyepiece:WF10X(Optional:WF16X);

Total Magnification:40X-1600X

Sliding Binocular Head,Interpupillary Distance:55-75mm;45°Inclined,360°

RotatableDouble Layers Mechanical Stage:140X155mm,Moving Range X-Y:70X50mm ϕ 2- ϕ 30mm Iris Diaphragm And ϕ 32

Filter:ABBE

Condenser:N.A.1.25

Coaxial Coarse And Fine Focusing Adjustable Mechanism:25mm,

Precision:0.002mm

Built-in illumination optional:Halogen Lamp 6V/20W,

Incandescent Lamp 15W,LED Lamp 1W,LEDlamp with Rechargeable Device

Net Weight/Set:6.4kg;

Package:1 Set/Wood Cabinet/Carton 38X34X50cm;

Gross Weight/Carton:11kg

