



AELAB BOD/COD Tester

Model:COD-100,



COD-100 Benchtop COD Analyzer



FEATURES

- COD analysis in accordance with EPA 410.4 and ISO 15705:2002 standards.
- LCD display screen.
- 2 LED light sources for different COD testing ranges.
- 2-point calibration.
- Data capacity of up to 200 sets GLP-compliant.
- Data stored can be transferred into the computer by USB communication interface.

COD-100R Reactor



Measuring slot
Convenient flip design



Meter includes
5pcs turbidity bottles

FEATURES

- 21 samples at a time.
- Settable heating temperature between 100 ~ 160 °C .
- Settable heating time between 0 ~ 120 min.
- Transparent plastic protection shield is equipped to protect the operator.

Model		COD-100
Measuring Unit	Range	(0 ~ 150)mg/L,(0 ~ 1500)mg/L
	Accuracy	±8%
	Repeatability	3%
	Fluctuation	6mg/L (20min)
General	Power Supply	AC Adapter, 200-240 VAC input
	Dimension (mm)& Weight (kg)	COD-100: 315×225×100, 2; COD-100R: 330×228×132, 4