

## Jar Tester & Automobile Exhaust Analyze

### Features:

- ◆ LCD display
- ◆ Synchro run, independent run



### Technical Specifications:

Model	ZR2	ZR4	ZR6
Beakers	2PCS	4pcs	6pcs
Acrylic Beaker	Round,1.0L, Φ11cm, H18cm Round,0.6L, Φ9.5cm, H13cm	Round, 0.6L, <P95*130mm	Round, 1.0L <P110*180mm Square optional, 1.5L, 110*110*170mm
Speed Range	10~1000rpm, stepless speed regulation		
Speed Accuracy	±0.5%		
Speed Gradient (G value)	10~1000/s. Just for round beaker Each step: 0~99min, 10 steps each program		
Running Time			
Measurement Temp. Range	0-50 °C		
Measurement Temp. Accuracy	±1 °C		
Consumption	40W	80W	180W
Power Supply	AC110/220V±10%, 50/60HZ		
External Size(W*D*H)	54×38×16mm	550*390*210mm	1080*350*400mm
Package Size( W*D*H)	56×47×25mm	610*450*270mm	950*480*530mm
Gross Weight	8.5KG	11kg	33kg