



CEB Micro Balance

High Capacity Balance

Model(型号)	CEB3C	CEB6C
Capacity(量程)	31g/3g	31g/6g
Readability(可读性)	0.01mg/0.001mg	0.01mg/0.001mg
Repeatability(重复性)	±0.02mg/±0.003mg	±0.02mg/±0.003mg
Linearity(线性)	±0.02mg/±0.010mg	±0.02mg/±0.010mg
Sensitivity drift(灵敏度)	±3 ppm/°C(10 to 30 °C/50 to 86 °F)	±3 ppm/°C(10 to 30 °C/50 to 86 °F)
Stabilization Time(稳定时间)	≤10s	≤10s
Accuracy Class(天平等级)	I	I
Min Weight(最小称量)	≥0.1mg	≥0.1mg
Operating environment(操作环境)	5 to 40 °C (41 to 104 °F), 85%RH or less	5 to 40 °C (41 to 104 °F), 85%RH or less
Product Size(Display unit)(产品尺寸:显示部分)	262*202*60mm	262*202*60mm
"Product Size(Weighing unit + breeze break unit)(产品尺寸:称重部分)"	325*138*107mm	325*138*107mm
"Net Weight(Display unit+Weighing unit + breeze break unit)(产品净重: 整体)"	4680g	4680g
Displayer(显示屏)	7" LCD Color Touch Screen	7" LCD Color Touch Screen
Unit(称重单位)	g, mg, oz, ct/customized	g, mg, oz, ct/customized
PC/Printer Function(电脑/打印机连接)	Yes	Yes
Communication interface(连接接口)	RS232, USB-b	RS232, USB-b
Plate Size(秤盘尺寸)	φ40mm	φ40mm
Plate Material(秤盘材质)	304SS	304SS
Volt(工作电压)	AC110-220V	AC110-220V
Frequency(工作频率)	50-60Hz	50-60Hz
Power(功率)	10W	10W

Electromagnetic Loadcell