

## BCS-14/35 Bacti-cinerator Sterilizer

Bacti-cinerator sterilizer is used for incineration and disinfection of small tools and bottles such as inoculation rings, inoculation needles, lab tweezers and triangle and tube mouth. It is easy to operate and widely used in the environment of biological safety cabinets, purification worktables, suction fans and mobile vehicles, etc.

### Features

- ◆ Replace alcohol lamp, safe and convenient
- ◆ Center temperature up to  $800^{\circ}\text{C} \pm 50^{\circ}\text{C}$ , rapid sterilization
- ◆ Small and light, easy to use
- ◆ Suitable for anaerobic chamber
- ◆ Ash organic material to prevent infectious splash and cross contamination
- ◆ Automatic shutdown after 3H, extending service life



### optional accessory



For BCS-14



For BCS-35



Fixed Angle



Max. Temp  $\geq 800^{\circ}\text{C} \pm 50^{\circ}\text{C}$



Rapid sterilize



Automatic shutdown

### Specifications

Model	BCS-14	BCS-35	Model	BCS-14	BCS-35
Central High Temperature	$800^{\circ}\text{C} \pm 50^{\circ}\text{C}$		Heater Head Length	150mm	100mm
Standby Temperature	$400^{\circ}\text{C} \pm 50^{\circ}\text{C}$		Power	200W	400W
Heater Head Angle	$21^{\circ}$		Weight	1.3kg	1.4kg
Max. Sample Diameter	$\phi 14\text{mm}$	$\phi 35\text{mm}$	Dimension (W×D×H)	187x112x200 (mm)	170x112x212 (mm)