

# TEXTILE BALANCE

## JA203SD

Count Balance

- ☆ Aluminum shell/Large space wind proof cover
- ☆ LCD(White back light with black font)
- ☆ Tare function
- ☆ Available in English and Chinese
- ☆ Overload alarm/Level indicator
- ☆ Interface RS232/Printer
- ☆ Date/Time/Number of samples/ Yarn blend codeset
- ☆ Automatically store 100 data/Delete and continuation
- ☆ Unit conversion:OZ/Tex/Nm/Ne/D/g/m<sup>2</sup>/oz/yd<sup>2</sup>
- ☆ Statistics/Print:Each weighing value, Maximum value(Max),Minimum value(Min), Average (AV),Uneven rate (CV)



JA203SD


### JA203SD specifications

Model	Capacity	Readability	Min Weighing	Repeatability	Linearity	Pan Size
JA203SD	210g	0.001g	0.004g	±0.002g	±0.003g	Φ90mm
/	Temperature	Interface	Draft Shield Size	Packing Size	Cal.Weight	Gross Weight
/	5-25°C	RS232	175*165*230mm	490*335*485mm	200g	9kg

## Y871

Daniel Balance

- ☆ LCD(White back light with black font)
- ☆ AC and DC power supply(AA\*3)
- ☆ Tare function/Counting/Unit conversion (g,OZ,Tex,Nm,Ne,D,g/m<sup>2</sup>,oz/yd<sup>2</sup>)
- ☆ Overload alarm/Level indicator
- ☆ Option: Dual display/Interface/Printer



XY-2CS

### Y871 specifications

Model	Capacity	Readability	Min Weighing	Pan Size	Temperatur	Cal. Weight	Packing size	Gross Weight
XY300-2CS	310g	0.01g	0.04g	Φ130mm	5-25°C	200g	310*230*155mm	2.5kg
XY600-2CS	610g	0.01g	0.04g	Φ130mm	5-25°C	200g	310*230*155mm	2.5kg