



DENSITY BALANCE

FA2204T SERIES

Density Balance

- ☆ Internal auto calibration/Electromagnetic force sensor
- ☆ One key/Time set internal calibration
- ☆ Aluminum shell/Large space wind proof cover
- ☆ Tare function/Counting/Percentage
- ☆ Under weighing/Unit conversion(g/mg/ct/oz...)
- ☆ Sensitivity/Speed set/Overload alarm/Fault alarm/Level indicator



FA2204T

Product specifications

Model	FA2204TS(solid)	FA2204TL(liquid)
Capacity	220g	
WR	0.0001g	
DR	0.0001g/cm ³	
Test Capacity	0.0001-99.9999g/cm ³	
Pan Size-01	60*152mm	Φ50mm
Pan Size-02	70*110mm	90*90*25mm
Tank Size	139*100*78mm WaterTank	50ml Beaker
Op-temp	15-25°C	
Display	Touch panel	
Cal. weight	Automatic internal calibration (default off) with 200g weight	
Material	Aluminium shell	
Outlook size	365*210*355 mm	
Packing Size	465*315*660mm	
Weight	11 kg	
Density Kit	Solid kit	Liquid Kit
Battery	No	

