


DENSITY BALANCE

JA-M SERIES

Density Balance



- ☆ Color plastic shell/Aluminum shell
- ☆ LCD/Tare function
- ☆ Overload alarm/Level indicator
- ☆ Interface:RS232(Standard)
- ☆ Option:RS485/Printer/Wind shield
- ☆ Solid/Liquid/Dust/Absorbent
- ☆ Density frame(Solid/Liquid/Powder)
- ☆ Direct display sample density
- ☆ Density reference set

Product specifications

| Mode Series | JA-M-S Solid JA-M-G Gold JA-M-A Absorbent | | | JA-M-L Liquid | | | JA-M-P Powder | | |
|---------------|---|--------------|--------------|-----------------|----------|----------|------------------------------|----------|----------|
| | JA203M-S/G/A | JA303M-S/G/A | JA602M-S/G/A | JA203M-L | JA303M-L | JA602M-L | JA203M-P | JA303M-P | JA602M-P |
| Model | JA203M-S/G/A | JA303M-S/G/A | JA602M-S/G/A | JA203M-L | JA303M-L | JA602M-L | JA203M-P | JA303M-P | JA602M-P |
| Capacity | 210g | 310g | 610g | 210g | 310g | 610g | 210g | 310g | 610g |
| WR | 0.001g | 0.005g | 0.01g | 0.001g | 0.005g | 0.01g | 0.001g | 0.005g | 0.01g |
| DR | 0.001g/cm ³ | | | | | | | | |
| Test Capacity | 0.001-99.999g/cm ³ | | | | | | | | |
| Pan Size-01 | 60*152mm | | | Φ50mm Upper pan | | | Φ90mm Weighing pan | | |
| Pan Size-02 | 70*110mm | | | 90*90*25mm | | | ----- | | |
| Tank Size | 139*100*78mm Water tank | | | 50ml Beaker | | | 50ml Specific gravity bottle | | |
| Op-temp | 5-35°C | | | | | | | | |
| Key name | Power,F,ENTER,CAL, 0 | | | | | | | | |
| Cal-weight | 200g | | | | | | | | |
| Material | Aluminium shell | | | | | | | | |
| Outlook size | 365*215*250mm | | | | | | | | |
| Packing Size | 410*270*165mm | | | | | | | | |
| Weight | 4.5kg | | | | | | | | |
| Density Kit | Solid kit | | | Liquid kit | | | Powder kit | | |
| Glass Draft | Option | | | | | | | | |
| Battery | No | | | | | | | | |