


DENSITY BALANCE

XY-2CM SERIES

Density Balance



- ☆ Color plastic shell
- ☆ LCD/Tare function
- ☆ Overload alarm/Level indicator
- ☆ Option:Interface(RS232/RS485) /Printer
- ☆ Solid/Liquid
- ☆ Density frame(Solid/Liquid)
- ☆ Direct display sample density
- ☆ Density reference set



Product specifications

Mode Series	XY-2CM-S Solid	XY-2CM-L Liquid
Model	XY3002CM-S	XY3002CM-L
Capacity	310g	
WR	0.01g	
DR	0.001g/cm ³	
Test Capacity	0.001-99.999g/cm ³	
Pan Size-01	60*152mm	Φ50mm
Pan Size-02	70*110mm	90*90*25mm
Tank Size	139*100*78mm	50ml Beaker
Op-temp	5-35°C	
Key name	Power,F,ENTER,CAL, 0	
Cal.weight	200g	
Material	Plastic shell	
Outlook size	270*195*230mm	
Packing Size	345*325*195mm	
Weight	4.5kg	
Density Kit	Solid kit	Liquid kit
Battery	Yes	

Density component

