

# MOISTURE ANALYZER

## XY-MW-T SERIES


Halogen Moisture Analyzer

PATENT NUMBER

201430389248.2



- ☆ 5 Inch touch screen/CN and EN switch
- ☆ Highlight halogen lamp/HBM sensor
- ☆ Aluminum shell/Stainless steel chamber
- ☆ Real-time dynamic curve/Date and Time
- ☆ Temperature/Time set
- ☆ The percentage of moisture content/dry residual
- ☆ Store historical set
- ☆ Option:Aluminum box/Printer/Grinder/USB





XY-MW-T


### Product application



MOISTURE MEASUREMENT REPORT	
MACHINE TYPE:	XY-100MW-T
HEATING METHOD	standard
STOP MODE	automatic
HEATING TEMPERATURE	130°C
HEATING TIME	04' 05"
WET WEIGHT	3.210g
DRY WEIGHT	2.815g
MOISTURE CONTENT	12.21%M



# MOISTURE ANALYZER



## Product specifications

Model	XY-100MW-T
Capacity	110g
Readability	0.001g
Operating Temp	5°C-35°C
Display	5 Inch touch screen
Calibration	External calibration
Interface	RS232(Standard) USB(Optional)
Pan Size	Φ90mm
Packing Size	485*340*350mm
Gross Weight	6.8kg
Heating Unit	Halogen lamp-2
Temp.Sensor	PT100
Temp.Range	40°C-199°C
Temp.Step	1°C
Moisture Range	0.00%-100.00%
Moisture Readability	0.01%
Dry Range	100.00%-0.00%
Dry Readability	0.01%
Temp Setting	40°C-199°C By 1°C Step
Time Setting	1-99min By 10s
Number Of Storage	20
Number Of History	20
Heater power	220V ±15% 50Hz / 110V±15% 60Hz / 400W
Balance power	Input:220V 50Hz / 110V 60Hz Output:9V