

MOISTURE ANALYZER

XY-MW SERIES


Halogen Moisture Analyzer

PATENT NUMBER

201430389248.2



- ☆ Halogen lamp/HBM sensor
- ☆ Aluminum shell/Stainless steel chamber
- ☆ LCD(White back light with black font)
- ☆ Temperature/Time set
- ☆ The percentage of moisture content/Dry residual
- ☆ Store historical set
- ☆ Option:Aluminum box/Printer/Grinder/USB



XY-MW

XY-MW-A SERIES

Halogen Moisture Analyzer



- ☆ Halogen lamp/HBM sensor
- ☆ Aluminum shell/Stainless steel chamber
- ☆ High definition LCD
- ☆ Highlight halogen lamp
- ☆ Temperature/Time set
- ☆ The percentage of moisture content/Dry residual
- ☆ Store historical set
- ☆ Option:Aluminum box/Printer/Grinder



MOISTURE ANALYZER



■ Product specifications

Model	XY-110MW	XY-105MW	XY-102MW	XY-100MW	XY-100MW-A
Capacity	110g				
Readability	0.01g	0.005g	0.002g	0.001g	
Operating Temp	5°C-35°C				
Display	Backlit LCD			High definition backlit LCD	
Calibration	External calibration				
Interface	RS232(Standard) USB(Optional)				
Pan Size	Φ90mm				
Packing Size	485*340*350mm				
Gross Weight	6.8kg				
Heating Unit	Halogen lamp-1			Halogen lamp-2	
Temp. Sensor	PT100				
Temp. Range	40°C-199°C				
Temp. Step	1°C				
Moisture Range	0.0%-100.0%	0.00%-100.00%			
Moisture Readability	0.2%	0.10%	0.04%	0.01%	
Dry Range	100.00%-0.00%				
Dry Readability	0.2%	0.10%	0.04%	0.01%	
Temp Setting	40°C-199°C By 1°C Step				
Time Setting	1-99min By 10s				
Number Of Storage	15				
Heate power	220V ±15% 50Hz / 110V±15% 60Hz / 400W				
Balance power	Input:220V 50Hz / 110V 60Hz Output:9V				
Chamber Height	22mm				