


HYDROSTATIC BALANCE

XY-J SERIES


Hydrostatic Balance



- ☆ For highway and textile industry
- ☆ Metal base
- ☆ Large size stainless steel plates
- ☆ LED(Green)
- ☆ Tare function/Counting
- ☆ Overload alarm
- ☆ Option:Interface/Printer



Product application



Product specifications

| Model | XY5000J | XY6000J |
|---------------|---------------|---------|
| Capacity | 5100g | 6100g |
| Readability | 0.1g | |
| Min Weighing | 0.4g | |
| Repeatability | ±0.2g | |
| Linearity | ±0.3g | |
| Pan Size | 160*160mm | |
| Temperature | 5-35°C | |
| Packing Size | 485*340*350mm | |
| Gross Weight | 6kg | |