

MOISTURE ANALYZER

XY-MW-IR SERIES

Infrared Moisture Analyzer

- ☆ Infrared lamp/HBM sensor
- ☆ Aluminum shell/Stainless steel chamber
- ☆ High definition LCD
- ☆ Temperature/Time set
- ☆ The percentage of moisture content/
Dry residual
- ☆ Store historical set
- ☆ Option:Aluminum box/Printer/Grinder



XY-MW-IR

Product specifications

Model	XY-100MW-IR
Capacity	110g
Readability	0.001g
Operating Temp	5°C-35°C
Display	High definition backlit LCD
Calibration	External calibration
Interface	RS232
Pan Size	Φ90mm
Packing Size	485*340*350mm
Gross Weight	7.0kg
Heating Unit	Infrared lamp
Temp.Sensor	PT100
Temp.Range	40°C-160°C
Temp.Step	1°C
Moisture Range	0.00%-100.00%
Moisture Readability	0.01%
Dry Range	100.00%-0.00%
Dry Readability	0.01%
Temp Setting	40°C-160°C By 1°C Step
Time Setting	1-99min By 10s
Number Of Storage	15
Number Of History	5
Heater power	220V ±15% 50Hz / 110V±15% 60Hz / 400W
Balance power	Input:220V 50Hz / 110V 60Hz Output:9V