

MOISTURE ANALYZER

XY-MS SERIES
Intelligent Moisture Analyze

DESIGN PATENT
2017304219332


UTILITY MODEL
2017211905798




High brightness halogen lamp

XY-MS

Product application



MOISTURE MEASUREMENT REPORT	
MACHINE TYPE:	XY-100MS(PLUS)
HEATING METHOD	standard
STOP MODE	automatic
HEATING TEMPERATURE	130°C
HEATING TIME	04' 05"
WET WEIGHT	3.210g
DRY WEIGHT	2.815g
MOISTURE CONTENT	12.21%M



MOISTURE ANALYZER

Product features

L*W*H(415*215*195 mm)

- ☆ 7 Inch touch screen
- ☆ Automatic door/Transparent window
- ☆ Highlight halogen lamp
- ☆ Polymer materials/Heat resistant molding
- ☆ Dual motor drive/Double fan design
- ☆ Temp.control/Stainless steel chamber
- ☆ Real-time dynamic curve/Date and Time
- ☆ Temperature/Time set
- ☆ Store historical set
- ☆ The percentage of moisture content
- ☆ Connect PC and printer
- ☆ Option:Printer/USB



Product specifications

Model	XY-100MS	XY-100MS(plus)	XY-200MS	XY-200MS(plus)
Capacity	110g		210g	
Readability	0.001g			
Operating Temp	5°C-35°C			
Pan Size	Φ100mm			
Heating Source	Halogen lamp			
Temp. Sensor	PT100			
Temp. Range	40°C-230°C			
Temp. Step	1°C			
Moisture Range	0.00%-100.00%			
Moisture Readability	0.01%			
Dry Range	100.00%-0.00%			
Dry Readability	0.01%			
Temp Setting	40°C-230°C By 1°C Step			
Time Setting	1-99min By 1s			
Number Of Storage	20			
Number Of History	20			
Heating Unit	220V ±15% 50Hz / 110V±15% 60Hz/400W			
Balance power	Input:220V 50Hz / 110V 60Hz Output:9V			
Machine power	Input: 220V 50Hz / 110V 60Hz Output:24V Dual-drive			
Host Dimension	415*215*195mm			
Indicator Dimension	200*135*80mm			
Operation Height	275mm			
The Chamber Height	35mm			
Display	7 Inch touch panel			
Interface	RS232(Standard)		RJ45/USB (optional)	
Communication Mode	Cable	Wireless	Cable	Wireless
Rechargeable Battery	No	10,000mAh	No	10,000mAh