

ANALYTICAL BALANCE

FA-A(I) SERIES

Internal Calibration Analytical Balance



Ioniser (Option)

- ☆ Internal calibration
- ☆ LCD(White back light with black font)
- ☆ One key/Time trigger internal calibration
- ☆ Electromagnetic force sensor
- ☆ Glass draft shield with 3 sliding doors for easy access to the items being weighed
- ☆ Tare function/Counting/Percentage
- ☆ Unit conversion(g/mg/ct/oz...)
- ☆ Sensitivity/Speed set
- ☆ Overload alarm/Fault alarm/Level indicator



FA (AI)

Product specifications

Model	FA314AI	FA414AI	FA514AI	FA614AI	FA1114AI
Capacity	310g	410g	510g	610g	1100g
Readability	0.1mg				
Min weighing	0.4mg				
Stable Time	≤4s				
Display	LCD(white back light with black font)				
Pan size	Φ80mm				
Operation temp	5-35°C				
Repeat ability	0.1mg if weight ≤ 220g	0.3mg from 220g to 310g	0.35mg from 310g to 610g		
Linearity	± 0.5mg				
Draft shield size	162mm × 171mm × 225mm				
Cal.weight	Internal calibration				
Interface	RS232(Standard) RS485/Printer(Option)				
Dimension	345 × 215 × 345mm				
N . weight	7 kg				