




ANALYTICAL BALANCE

FA-T SERIES

Internal Calibration Touch Screen Analytical Balance

L*W*H(365*210*355mm)



Product features

- ☆ Internal auto calibration
- ☆ 5 Inch touch screen
- ☆ One key/Time trigger internal calibration
- ☆ Electromagnetic force sensor
- ☆ Aluminum shell/Large space wind proof cover
- ☆ Tare function/Counting/Density/Percentage
- ☆ Unit conversion(g/mg/ct/oz...)
- ☆ Under weighing
- ☆ Sensitivity/Speed set
- ☆ Overload alarm/Fault alarm/Level indicator
- ☆ Option:Density kit (Solid/Liquid)



Product application



Gold jewelry



High precision data matching



Scientific experiment




Medicine

ANALYTICAL BALANCE

FA-N SERIES

Internal Calibration Analytical Balance

L*W*H(365*210*355mm)



- ☆ Internal auto calibration
- ☆ LCD(White back light with black font)
- ☆ One key/Time set internal calibration
- ☆ Electromagnetic force sensor
- ☆ Aluminum shell/Large space wind proof cover
- ☆ Tare function/Counting/Percentage
- ☆ Unit conversion(g/mg/ct/oz...)
- ☆ Under weighing
- ☆ Sensitivity/Speed set
- ☆ Overload alarm/Fault alarm/Level indicator



FA-N

Product specifications

Model	FA1004T	FA1104T	FA1204T	FA2004T	FA2104T	FA2204T
Model	FA1004N	FA1104N	FA1204N	FA2004N	FA2104N	FA2204N
Capacity	100g	110g	120g	200g	210g	220g
Readability	0.0001g					
Min Weighing	0.0004g					
Repeatability	±0.0002g					
Linearity	±0.0003g					
Stable Time	≤3S					
Pan Size	Φ80mm					
Appearance Size	365*210*355mm					
Draft Shield Size	175*165*230mm (Inner diameter)					
Communication	RS232 (Standard) USB (Option)					
Temperature	15-35°C					
Calibration	Internal calibration					
Cal.Weight	100g			200g		
Option	Printer					
Packing Size	Inside:440*295*430mm Outside:490*335*485mm					
Gross Weight	11kg					