

MOISTURE ANALYZER

XY-MW-R SERIES

Moisture Regain Analyzer

PATENT NUMBER
201430389248.2


- ☆ Moisture regain of textile
- ☆ Test cotton, yarn, cloth, cashmere, silk, chemical fibrils...
- ☆ Stainless steel net-cover
- ☆ Highlight halogen lamp
- ☆ HBM sensor
- ☆ Aluminum shell
- ☆ Stainless steel chamber
- ☆ LCD(White back light with black font)
- ☆ Temperature/Time set
- ☆ The percentage of moisture regain
- ☆ Store historical set
- ☆ Option:Aluminum box/Printer



XY-MW-R

Product specifications

1. Stainless steel upper mesh cover
2. Stainless steel lower mesh cover
3. Stainless steel plate



Model	XY-105MW-R	XY-100MW-R
Capacity	110g	
Readability	0.005g	0.001g
Operating Temp	5°C-35°C	
Display	LCD	
Calibration	External calibration	
Interface	RS232	
Net Cover Size	Φ120mm	
Packing Size	485*340*350mm	
Gross Weight	7.0kg	
Heating Unit	Halogen lamp	
Temp.Sensor	PT100	
Temp.Range	40°C-199°C	
Temp.Step	1°C	
Regain Range	0.00%-+∞	
Moisture Readability	0.10%	0.01%
Temp Setting	40°C-199°C By 1°C Step	
Time Setting	1-99min By 10s	
Number Of Storage	15	
Heater power	220V ±15% 50Hz /110V±15% 60Hz/400W	
Balance power	Input:220V 50Hz/110V 60Hz Output:9V	