

PRECISION BALANCE

JA-C SERIES

Precision Balance

L*W*H(345*205*105mm)




Product function

- ☆ Aluminum shell
- ☆ Unit conversion(g/ct/oz...)
- ☆ LED (Green) /LCD
- ☆ Tare function/Counting
- ☆ Overload alarm/Level indicator
- ☆ Option:Printer



Product display



Product specifications

Model	JA2000C	JA3000C	JA4000C	JA5000C
Capacity	2000g	3000g	4000g	5000g
Readability	0.01g			
Repeatability	±0.02g			
Linearity	±0.03g			
Pan Size	180*180mm			
Stable Time	≤3S			
Interface	RS232(Standard)			
Dimensions	345*205*105mm			