

PRECISION BALANCE

JA-H SERIES

Precision Balance

L*W*H(365*210*355mm)



Product function

- ☆ Aluminum shell / Large space wind proof cover
- ☆ LED(Green)/LCD
- ☆ Tare function/Counting
- ☆ Overload alarm/Level indicator
- ☆ Option:Under weighing/Printer



JA-H(LED)



JA-H(LCD)

Product specifications

Model	JA103H	JA203H	JA303H
Capacity	110g	210g	310g
Readability	0.001g		
Min Weighing	0.004g		
Repeatability	±0.002g		
Linearity	±0.003g		
Stable Time	≤3S		
Pan Size	Φ90mm		
Appearance Size	365*210*355mm		
Draft Shield Size	175*165*230mm (Inner diameter)		
Interface	RS232(Standard)		
Calibration	External calibration		
Cal.Weight	100g	200g	
Option	Printer/Under weighing		
Packing Size	Inside:440*295*430mm Outside:490*335*485mm		
Gross Weight	7.5kg		