

PRECISION BALANCE

JA-P SERIES

Precision Balance

L*W*H(270*200*280mm)



Product function

- ☆ LCD/AC and DC power supply(AA*3)
- ☆ Tare function/Counting/Unit conversion(g/ct/oz...)
- ☆ Min weighing set
- ☆ Overload alarm/Level indicator
- ☆ Option: Dual display/Interface/Printer/Dust cover



JA-P

*Optional extra display

Product specifications

Model	JA103P	JA203P	JA303P	JA103PL
Capacity	110g	210g	310g	110g/610g
Readability	0.001g			0.001g/0.01g
Min Weighing	0.001g-0.004g			1d-4d
Repeatability	±0.002g			±2d
Linearity	±0.003g			±3d
Stable Time	≤3S			
Pan Size	Φ90mm			Φ110mm
Appearance Size	270*200*280mm			
Draft Shield Size	170*165*180mm (Inner diameter)			
Calibration	External calibration			
Cal.Weight	100g	200g		
Packing Size	360*265*350mm / 560*380*730mm(4 units per box)			
Gross Weight	4kg / 16kg(4 units per box)			