

PRECISION BALANCE


XY-1B/1BF SERIES

Precision Balance

L*W*H(190*230*70mm)



- ☆ Fashional appearance
- ☆ Large size stainless steel pan(For XY-1BF)
- ☆ LCD(White back light with black font)
- ☆ AC and DC power supply(AA*3)
- ☆ Tare function/Counting/Unit conversion(g/ct/oz...)
- ☆ Min weighing set
- ☆ Overload alarm/Level indicator
- ☆ Option: Dual display/Interface/Printer/Dust cover Calibration weight



Product specifications

Model	Capacity	Readability	Repeatability	Linearity	Pan Size	Stable Time
XY300-1B	310g	0.1g	±0.2g	±0.3g	Φ130mm	≤2S
XY600-1B	610g	0.1g	±0.2g	±0.3g	Φ130mm	≤2S
XY1000-1B	1100g	0.1g	±0.2g	±0.3g	Φ130mm	≤2S
XY2000-1B	2100g	0.1g	±0.2g	±0.3g	Φ130mm	≤2S
XY1000-1BF	1100g	0.1g	±0.2g	±0.3g	160x160mm	≤2S
XY2000-1BF	2100g	0.1g	±0.2g	±0.3g	160x160mm	≤2S
XY3000-1BF	3100g	0.1g	±0.2g	±0.3g	160x160mm	≤2S
XY5000-1BF	5100g	0.1g	±0.2g	±0.3g	160x160mm	≤2S
XY6000-1BF	6100g	0.1g	±0.2g	±0.3g	160x160mm	≤2S
XY10K-1BF	10000g	0.1g	±0.2g	±0.3g	160x160mm	≤2S

Product application

