

PRECISION BALANCE

XY-2C SERIES

Precision Balance

L*W*H(190*230*70mm)



- ☆ LCD/AC and DC power supply(AA*3)
 - ☆ Tare function/Counting/Unit conversion(g/ct/oz...)
 - ☆ Min weighing set
 - ☆ Overload alarm/Level indicator
 - ☆ Option: Dual display/Interface/Printer/ Windproof/Dust cover
- Note: capacity more than 1kg or accuracy less than 0.01g is not equipped with calibration weight



■ Product specifications

Model	Capacity	Readability	Repeatability	Linearity	Pan Size	Stable Time
XY100-2C	110g	0.01g	±0.02g	±0.03g	Φ130mm	≤2S
XY200-2C	210g	0.01g	±0.02g	±0.03g	Φ130mm	≤2S
XY300-2C	310g	0.01g	±0.02g	±0.03g	Φ130mm	≤2S
XY500-2C	510g	0.01g	±0.02g	±0.03g	Φ130mm	≤2S
XY600-2C	610g	0.01g	±0.02g	±0.03g	Φ130mm	≤2S
XY1000-2C	1100g	0.01g	±0.02g	±0.03g	Φ130mm	≤2S
XY2000-2C	2100g	0.01g	±0.02g	±0.03g	Φ130mm	≤2S
XY3000-2C	3100g	0.01g	±0.02g	±0.03g	Φ130mm	≤2S