

PRECISION BALANCE

XY-EB SERIES

Precision Balance



- ☆ LCD display
- ☆ Stainless steel weighing pan (Both side can be used)
- ☆ USB charging port (No adapter)
- ☆ 2 x AAA dry battery
- ☆ Units conversion: g/oz/ct/tl



XY-EB

Product application



Product specifications

| Model | XY3KB1EB | XY6KB5EB | XY10KA1EB |
|-----------------|-----------------------------|----------|-----------|
| Capacity | 3000g | 6000g | 10000g |
| Readability | 0.1g | 0.5g | 1g |
| Repeatability | 0.2g | 1.0g | 2g |
| Linearity | 0.3g | 1.5g | 3g |
| Min Weighing | 0.4g | 2.0g | 4g |
| Stable Time | ≤3s | | |
| Pan Size | 157*157mm | | |
| Temperature | 5-35°C | | |
| Appearance Size | 160*225mm | | |
| Net weight | 0.6kg per unit/15kg per box | | |
| Gross weight | 0.8kg per unit/17kg per box | | |