

# TEXTILE BALANCE

## YG101

Raw Silk Size Tester



- ☆ Aluminum shell/Large space wind proof cover
- ☆ LCD (White back light with black font)
- ☆ Tare function
- ☆ Language: Chinese / English
- ☆ Overload alarm/Level indicator
- ☆ Standard: Interface RS232/Printer
- ☆ Date/Time/Number of samples set
- ☆ Automatically store 200 data/ Delete and continuation
- ☆ Inquiry/Print: Testing date,time,checker, number of machinery,the distribution of raw silk size,hank number,sum,Average fineness, Denier deviation,maximum diviation



YG101

### YG101 specifications


Model	Capacity	Readability	Min Weighing	Repeatability	Linearity	Pan Size
YG101	110g(99.5D)	0.001g(0.5D)	0.004g	±0.002g	±0.003g	Φ90mm
/	Temperature	Interface	Draft Shield Size	Packing Size	Cal.Weight	Gross Weight
/	5-25°C	RS232	175*165*230mm	490*335*485mm	100g	9kg

## Y871

Daniel Balance



- ☆ LED(Green)
- ☆ Tare function/Counting
- ☆ Overload alarm
- ☆ Option: Under weighing/Interface/Printer
- ☆ Number of samples set
- ☆ Display:  
G (g) / Dan Neal value (D) (size value)



Y871

### Y871 specifications

Model	Capacity	Readability	Min Weighing	Pan Size	Temperatur	Cal. Weight	Packing Size	Gross Weight
Y871	110g(99.5D)	0.001g(0.5D)	0.004g	Φ130mm	5-25°C	100g	310*230*155mm	3kg