

SEMI-MICRO ANALYTICAL BALANCE

FA(SEM) SERIES

Semi-micro Analytical Balance



Ioniser (Option)

- ☆ Recipe database: 99 datas can be stored
- ☆ Automatic internal calibration
- ☆ Tare/counting/Percentage test/Animal weighing
- ☆ Baud rate is adjustable/ Density function
- ☆ Unit conversion (g/mg/ct/oz...)
- ☆ Multi-Language: English/ German/French/ Italian/Spanish/Portuguese
- ☆ Keep weight values in GLP/ISO format, Includes date, time, and serial number
- ☆ Plastic protection cover for keypad
- ☆ Overload alarm/Fault alarm/Level indicator



FA (SEM)

Product specifications

Model	FA1265SEM	FA2285SEM	FA22105SEM	FA125SEM-ION	FA165SEM-ION	FA205SEM-ION
Capacity	62g/120g	82g/220g	102g/220g	122g	162g	210g
Readability	0.01mg/0.1mg			0.01mg		
Min weighing	60mg		100mg	100mg		
Stable Time	≤6s					
Display	LCD(white back light with black font)					
Pan size	Φ80mm					
Operation temp	5-35°C 20-80% Humidity					
Repeat ability	0.03mg/0.1mg	0.05mg/0.1mg	0.05mg	0.05mg		
Linearity	±0.08mg/±0.2mg	±0.1mg/±0.2mg	±0.1mg	± 0.1mg		
Draft shield size	152mm × 160mm × 225mm					
Cal.weight	Internal calibration					
Option	Printer					
Dimension	345 × 215 × 345mm					
G.N . weight	8.25kg/5.8kg			8.3kg/5.85kg		
IONISER	Option			Standard		
Language	Italian, English, Portugues, Deutsch, France, Espanol					
Packing Size	510x400x620mm					
Power Supply	110-230V 50/60Hz Out DC 12V 1A					