

AE LAB Water Bath

- ◆ Microprocessor control,
- ◆ Digital display inner temperature
- ◆ Stainless steel interior
- ◆ The aluminum ring cover
- ◆ Suitable for thermostatic , baking,dewaxing, sterilization



HH-S4

Model	HH-S1	HH-S2	HH-S4	HH-S6	HH-S8
Capacity	3L	6L	10L	16L	22L
Power	300W	500W	1000W	1500W	2000W
Temperature fluctuation	±0.5℃	±0.5℃	±1℃	±1.5℃	±1.5℃
Temperature range	RT+5-99℃				
Temperature homogeneity	±1℃				
Voltage	110V 60Hz/220V 50Hz				
Working size(mm)	160x160x120	300x160x120	300x300x120	460x300x120	600x300x120
Dimension (mm)	230x230x240	370x230x240	370x370x240	530x370x240	670x370x240
Package (mm)	290x290x280	410x280x280	420x410x280	610x430x280	730x430x280
G.W.IN.W.	3.5/3kg	4.5/4kg	6/5kg	7.5/6.5kg	9/8kg

- ◆ Microprocessor control,
- ◆ Digital display inner temperature
- ◆ One-timer-molding cistem,
- ◆ The aluminum ring cover
- ◆ Stainless steel roof cover (optional)
- ◆ Suitable for thermostatic , baking, dewaxing,sterilization



WH-4

Model	WH-1	WH-2	WH-4	WH-6
Capacity	3L	6L	12L	16L
Power	200W	500W	1000W	1500W
Temperature fluctuation	±0.5℃	±0.5℃	±1℃	±1.5℃
Temperature range	RT+5-99℃			
Temperature homogeneity	±1℃			
Voltage	110V 60Hz/220v 50Hz			
Working size(mm)	150x140x120	300x150x120	320x320x120	460x300x120
Dimension (mm)	180x180x240	330x180x240	360x330x240	540x330x250
Package (mm)	230x220x300	380x230x280	410x400x280	590x400x290
G.W./N.W.	3/2.5kg	4/3kg	6/5kg	8/7kg

- ◆ Microprocessor control,
- ◆ Digital display setting and actual temperature
- ◆ Digital display running timer
- ◆ PID self-tuning , over temperature alarm
- ◆ Automatic alarm after operation
- ◆ Suitable for thermostatic , baking, dewaxing,sterilization


DK-S26

DK-8A
DK-600B

- ◆ Microprocessor control,
- ◆ Digital display setting and actual temperature and time
- ◆ PID self-tuning , over temperature alarm
- ◆ Automatic alarm after operation
- ◆ Good insulation performance with stainless steel roof cover
- ◆ With bakelite handle
- ◆ DK-600B with double layer cover.

Model	DK-S12	DK-S24	DK-S26	DK-S28
Capacity	6L	10L	16L	22L
Power	500W	1000W	1500W	2000W
Temperature fluctuation	±0.5℃	±0.5℃	±1℃	±1℃
Temperature range	RT+5 -99℃			
Temperature homogeneity	±1℃			
Timer	0-999min			
Voltage	110V 60Hz /220V 50Hz			
Working size(mm)	300x160x120	300x300x120	460x300x120	600x300x120
Dimension(mm)	470x230x240	470x370x240	630x370x240	770x370x240
Package(mm)	540x290x280	540x430x280	710x430x280	880x430x280
G.W./N.W.	5.5/5kg	7/6.5 kg	11/10 kg	13/12 kg

Model	DK-320	DK-420	DK-8A	DK-600B
Capacity	6L	9L	23L	20L
Power	500W	1000W	1200W	1500W
Temperature fluctuation	±0.5℃	±0.5℃	±1℃	±1℃
Temperature range	RT+5 -99℃			
Temperature homogeneity	±0.5℃			
Timer	0-999min			
Voltage	110V 60Hz /220V 50Hz			
Working size(mm)	300x160x120	400x180x130	450x300x170	580x280x130
Dimension(mm)	470x230x260	570x270x280	630x370x320	760x360x260
Package(mm)	530x300x310	630x320x320	670x430x360	810x420x320
G.W./N.W.	7/6kg	8/7kg	12/11kg	15/14kg

- ◆ Microprocessor control,
- ◆ Digital display inner temperature
- ◆ Good insulation performance with stainlesssteel roof cover
- ◆ With bakelite handle


HH-W420

WB-18S

- ◆ Microcomputer PID intelligent control and LCD backlight display temperature and time
- ◆ The upper cover is equipped with ring butterfly cover
- ◆ Stainless steel inner gallbladder tensile forming process,the inner wall is smooth, the product is more beautiful
- ◆ Over temperature protection, leakage current protection,low water alarm protection
- ◆ Prevent dry burning heating components, ffectively protect heating elements
- ◆ Optional stainless steel flat cover or roof cover

Model	HH-W420	HH-W600
Capacity	9L	20L
Power	1000W	1500W
Temperature range	RT+5 -99℃	RT+5 -99℃
Temperature fluctuation	±1℃	±1℃
Temperature homogeneity	±1℃	±1℃
Voltage	110V 60Hz /220V 50Hz	110V 60Hz /220V 50Hz
Working size(mm)	400x180x130	580x280x130
Dimension(mm)	480x260x280	660x360x280
Package(mm)	530x320x320	710x420x330
G.W./N.W.	7/6kg	11/10kg

Model	WB-5S	WB-12S	WB-18S
Volume	5L	12L	18L
Working Size (mm)	300x150x120	320x320x120	460x300x130
Power	500W	1000W	1500W
Temperature range	RT+5-99℃	RT+5-99℃	RT+5-99℃
Temperature accuracy	±0.3℃	±0.3℃	±0.3℃
Timer	0-9999Hour	0-9999Hour	0-9999Hour
Voltage	220V 50Hz	220V 50Hz	220V 50Hz
Dimension(mm)	350x200x270	370x370x270	530x320x290
Package(mm)	400x260x320	430x430x320	590x380x350
G.W./N.W.	6/5kg	10/9kg	13/12kg