

Electric Centrifugal Divider LXFY-2

Features

The Precision Divider has repeatedly demonstrated its ability to precisely divide larger samples into working size samples in a matter of seconds. It is ideal for smaller seeded crops. Dividing and mixing are accomplished by the centrifugal action of a motor driven revolving neoprene disc positioned under the hopper. A smooth inside chrome plated finish virtually eliminates any of the samples sticking to the inner surfaces. Both the hopper and cylinder top lift off to provide easy access to the interior of the dividing chamber.



Technical parameter

Hopper capacity	2 kg
Work error	< $\pm 0.5\%$ per kg
Size	33*33*49cm
Packing Size	63*40*37cm
Power	220V 50HZ
Packing list	One unit of machine and 3 units of sample holder