

LDO Forced Air Drying Oven

Features

- ★ Forced air, good temperature uniformity.
- ★ Microprocessor controller (with temperature correction and timing function), accurate and reliable.
- ★ Large LCD screen display.
- ★ High-quality stainless-steel chamber, removable shelf, easy-to-clean. Powder coating steel housing.
- ★ Equipped with leakage protection and a spare mechanical thermostat to ensure the product work normally even the main temperature control failed (for heating).
- ★ Anti-hot handle



Options:

2. Multi-segment programmable control
4. Built-in Printer
3. RS485 interface

Technical Parameters:

- 1.Temp. Fluctuation: $\pm 1^{\circ}\text{C}$.
- 2.Temp. Uniformity: $\pm 2.5^{\circ}\text{C}\%$.
- 3.Temp.Resolution: 0.1°C
- 4.Timer: 1~9999 mins



Specifications:

Model	Volume	Chamber Size (WxDxH) cm	Exterior Size (WxDxH) cm	Power Rating	Temp. Range	Shelf	Power Supply
LDO-9030A	31.5L	30x30x35	45x47x70	700W	RT+5~ 200°C	2/5	AC 220V 50HZ
LDO-9070A	72L	40x40x45	55x60x80	1000W		2/7	
LDO-9140A	136L	45x55x55	64x73x90	1600W		3/8	
LDO-9240A	225L	50x60x75	88x72x93	2000W		3/11	AC 380V 50HZ
LDO-9420A	429L	60x55x130	78x71x170	2600W		4/6	
LDO-9620A	624L	80x60x130	98x80x178	3600W	5/6	50HZ	
LDO-9036A	31.5L	30x30x35	45x47x70	950W	RT+5~ 300°C	2/5	AC 220V 50HZ
LDO-9076A	72L	40x40x45	55x60x80	1400W		2/7	
LDO-9146A	136L	45x55x55	64x73x90	2000W		3/8	AC380V 50HZ
LDO-9246A	225L	50x60x75	69x77x114	2600W		3/11	
LDO-9426A	429	60x55x130	78x71x170	3500W		4/6	
LDO-9626A	624	80x60x130	103x80x178	5700W		5/6	50HZ

★ LDO-9240A/ LDO-9246A and LDO-9620A/ LDO-9626A equipped with a FAN switch