

Hot Air Sterilizer

Features:

- ★ Forced air convection, good temperature uniformity.
- ★ Microprocessor controller (with temperature correction and timing function), accurate and reliable.
- ★ Large (2.2 inch) LCD screen display
- ★ High-quality stainless-steel chamber, removable shelf, easy-to-clean.
- ★ Equipped with leakage protection and a spare mechanical thermostat to ensure the product work normally even the main temperature control failed.
- ★ Designed with a transparent window made of reinforced tempered glass, easy observation.
- ★ Anti-hot handle



Options:

10. Multi-segment programmable control
11. Built-in Printer
12. RS485 interface

Specifications:

1. Temp. Fluctuation: $\pm 1^{\circ}\text{C}$
2. Temp. Uniformity: $\pm 2.5^{\circ}\text{C}\%$.
3. Temp. Resolution: 0.1°C
4. Timer: 1~9999mins



Technical Specifications:

Model	LAS-9023A	LAS-9053A	LAS-9073A	LAS-9123A	LAS-9203A
Volume	25L	54L	81L	138L	252L
Voltage	220V, 50HZ				
Temp. Range	RT+10~250 $^{\circ}\text{C}$				
Timing Range	1~9999min				
Operating Temp.	+5~40 $^{\circ}\text{C}$				
Power Rating	650W	1000W	1500W	1600W	2450W
Chamber Size (W×D×H)	34×30×32cm	42×33×35cm	45×40×45cm	50×50×55cm	60×60×70cm
Exterior Size ((W×D×H)	63×45×52cm	71×52×55cm	73×53×69cm	78×53×80cm	88×73×95cm
Grid Racks (Std/Max)	2/4	2/5	2/6	3/7	3/10