

Portable Hardness Tester

Features:

- *Measure wall and top thickness of various capsules.
- * High quality precision per centimeter measurement, high testing accuracy.

Technical Specifications:

Model	HD-1	HD-2	HD-3	HD-4
Thickness test range	0-5 mm	0-30 mm	0-5 mm	0-30 mm
Thickness measurement accuracy	±0.01mm			
Measurable capsule type	Capsule 00, 0, 1, 2, 3, 4, 5			
Dimensions	200*80*50 mm		160*60*40 mm	
Weight	0.9kg		0.3kg	




CLARIFYTESTER

ADVANTAGE:

- Flash eliminated. Illumination is adjustable.
- Timing and alarming function automatically. Test continuously.
- The three-primary-color fluorescent lamp is controlled by good lighting sensor, electronic inductor and shading.

Specification



Model	YB-2A	YB-3
Type	One side for one operator	Two sides for two operators
Lamp	Three-primary-color fluorescent lamp	
Illumination Range(lux)	1000—6000LX	
Time Range(s)	1~99	
Power Supply	AC110V/220V±10%, 50/60Hz	
Consumption	50W	
Packing Size(mm)	430*320*200mm	670*280*720mm
Gross Weight(kg)	6.5	11