

Soxhlet Extraction Fat Analyzer

Fat Analyzer tests fat content according to Soxhlet Extraction principle It features extraction, reagent recovery, cooling system and temperature controlling.

Features:

- ◆ Measuring range: oil content within the scope of 0.5%-60% of good, feed, oil and other fat products.
- ◆ Measuring sample quantity: 6 simultaneously.
- ◆ Heating time: within 10 min.
- ◆ Recovery system: automatic.
- ◆ Recovery rate: N80%.
- ◆ Temperature controlling range: room temperature -100°C . Power: AC220V/50HZ, 1000W.



Technical Specifications:

Model	SZF-06C
Shape size	670x250x600mm
Weight	23kg
Difference	Even heating and automatic recycling, automatic soaking and extraction, save working time.