



AELAB Water Bath

The water bath is a based heating equipment, widely used in enterprises, hospitals, colleges, research institutes and government laboratories.

Features:

- ◆ LCD screen shows data easy to understood and operate.
- ◆ PID temperature control accuracy.
- ◆ Over-temperature alarm, safe and reliable.
- ◆ The chamber is made of stainless-steel (304).
- ◆ Water bath has five covers which diameter size is 30mm, 50mm, 70mm, 95mm and 115mm.



Specifications:

Model	HWS-11	HWS-12	HWS-24	HWS-26	HWS-28
Hole	One	Two	Four	Six	Eight
Power Requirement	AC220V 50HZ				
Temperature Range	RT+5~100℃				
Temperature Stability	±0.5℃				
Display Resolution	0.1℃				
Timing Range	1~9999min				
Power Consumption	400W	500W	1000W	1500W	2000W
Interior Dimension (W×D×H mm)	155×155×110	305×155×110	305×305×110	455×305×110	600×300×110
External Dimension (W×D×H mm)	195×195×220	340×195×220	345×345×220	495×345×220	780×380×180
Volume	2.5L	4.9L	9.9L	14.8L	19.8L

Remark:

1. The technical parameters are tested at the conditions of ambient temperature 20℃, 50%RH, no load, no strong magnetic and no vibration.
2. The picture is for reference only, please refer to the actual goods.