

## Single-Channel Manual Pipettes



The Single-Channel Manual Pipette provides accurate liquid transfer with a lightweight, ergonomic design. Its corrosion-resistant material supports autoclaving at 121°C, and the elastic suction nozzle ensures leak-proof operation. One-button locking and smooth pipetting enhance usability and precision.

### Features:

- Swift & Precise
- Ergonomic Design
- Flexible tip
- Smart Lock
- Autoclavable
- Smooth Pipetting

### Specifications

Item Number	Nominal Volume ( $\mu$ l)	Step Volume ( $\mu$ l)
602021	0.5 - 10	0.01
602051	5 - 100	0.1
602061	10 - 300	1
602071	50 - 1000	1