

Halogen Moisture Analyzer

Features

- 5.5" LCD Touch screen.
- Aluminum alloy base and ABS shell.
- Timer dry / Auto dry function.
- RS232 interface, support to connect with printer and PC.



Specification

Model	AE-HM005	AE-HM105	AE-HM102	AE-HM101
Capacity	50g	100g	100g	100g
Readability	0.005g	0.005g	0.002g	0.001g
Moisture Range	0~100%	0~100%	0~100%	0~100%
Accuracy for Moisture measure	±0.5% (Sample≥5g)	±0.5% (Sample≥5g)	±0.5% (Sample≥5g)	±0.5% (Sample≥5g)
Moisture Readability	0.02%	0.02%	0.01%	0.01%
Repeatability	0.2%	0.2%	0.1%	0.0%
Interface	None	RS232 Interface	RS232 Interface	RS232 Interface
Scale Pan Size	Φ100mm	Φ100mm	Φ100mm	Φ100mm
Heating Way	High efficiency Halogen Lamp			
Heating Time	1~99 minutes interval of 1 second			
Time End Control mode	Auto Stop / Timed Stop			
Sample Drying Record	8 times	8 times	8 times	8 times
Temperature Range	60~160°C	60~160°C	60~160°C	60~160°C
Temp. Increase Step	±1°C	±1°C	±1°C	±1°C
External Size (WDH)	205x325x170mm	205x325x170mm	205x325x170mm	205x325x170mm
Shipping Size (WDH)	325x460x295mm	325x460x295mm	325x460x295mm	325x460x295mm
Net Weight	2.9kg	2.9kg	2.9kg	2.9kg
Gross Weight	4.2kg	4.2kg	4.2kg	4.2kg